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### IMPORTANT SAFETY INSTRUCTIONS

- Digital medical thermometer MT 1622 is intended to be used for oral, axillary or rectal measuring of human temperature. Field of application: in medical institutions and at home for individual use. A display of values with an accuracy of 0.01 °C. It can be recommended for family planning.
- The instrument may be used only for measuring body temperature!
- The minimum measurement time until the signalling tone [beep] must be maintained without exception!
- The thermometer contains small parts (battery, etc.) that can be swallowed by children. Therefore never leave the thermometer unattended to children.
- Protect the thermometer from hitting and striking!
- Avoid ambient temperatures above 60 °C. NEVER boil the thermometer!
- Use only the commercial disinfectants listed in the section «cleaning and disinfecting» to clean the thermometer. The thermometer must be intact when immersed in liquid disinfectant.
- The manufacturer recommends to verify the accuracy by an authorized laboratory every 2 years.

Batteries and electronic instruments must be disposed of in accordance with the locally applicable regulations, not with domestic waste.

### TURNING ON THE THERMOMETER

To turn on the thermometer, push the button over the display window; a short signalling tone [beep] signals «thermometer ON». A display test is performed. All signals should be displayed. Then at an ambient temperature of less than 32 °C, an «L» and a flashing «°C» appear at the upper right of the display. The thermometer is now ready for use.

### USING THE THERMOMETER

Choose the preferred measuring method. During a measurement, the current temperature is displayed continuously. The «°C» symbol flashes during this period of time. If the signalling tone sounds [beep-beep 10 times] and the «°C» is no longer flashing, this means that the measured increase in temperature is now less than 0.01 °C in 8 seconds and that the thermometer is ready for reading.

To prolong the life of the battery, turn it off after use by pressing the operating button briefly. Otherwise the thermometer will automatically turn off after about 10 minutes.

### STORAGE OF MEASURED VALUES

If you hold the button down for more than 3 seconds when turning the thermometer on, the maximum temperature stored automatically at the last measurement will be displayed. At the same time, a small «M» for memory will appear at the right of the display. After 3 more seconds, this value will disappear and the thermometer will switch to the normal measurement mode. Then the stored value is erased.

### BATTERY REPLACEMENT

When the «V»-symbol appears, battery is exhausted and needs replacing. Please remove the compartment lid of the thermometers and replace the exhausted battery by positioning the + at the top.

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This device complies with the requirements of the Medical Device Directive 93/42/EEC. Technical alterations reserved.

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### GUARANTEE

We grant you a 5 year guarantee after the date of purchase. Any damage caused by improper handling shall not be covered by the warranty. Battery and packaging are also excluded from the warranty. All other damage claims excluded. A warranty claim must be submitted with the purchase receipt. Please pack your defective thermometer well and send with sufficient postage to:

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Microlife AG, Espenstrasse 139, 9443 Widnau, Switzerland